

Product Datasheet

Eplerenone (orb1310957)

Description	Eplerenone
Storage	Store at -20°C for 12 months.
Note	For research use only
Formula	$C_{24}H_{30}O_6$
Purity	99.44%
MW	414.49
CAS Number	107724-20-9
SMILES	<chem>C[C@@]12[C@]34[C@]([C@]5([C@](C)(C[C@]3(O4)[H])[C@@]6(CC5)CC(=O)O6)[H])([C@H](C(OC)=O)CC1=CC(=O)CC2)[H]</chem>
Expiration Date	12 months from date of receipt.