


Product Datasheet

NSE (ENO2) antibody (orb1311066)

Description	Mouse monoclonal antibody to NSE (ENO2)
Species/Host	Mouse
Reactivity	Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 188-293 of human NSE (NP_001966) produced in E.coli.
Target	ENO2
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	47.3 kDa
Uniprot ID	P09104
NCBI	NM_001975
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...