


## Product Datasheet

### NSE (ENO2) antibody (orb1311236)

<b>Description</b>	Mouse monoclonal antibody to NSE (ENO2)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 188-293 of human NSE (NP_001966) produced in E.coli.
<b>Target</b>	ENO2
<b>Note</b>	For research use only
<b>Application notes</b>	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	47.3 kDa
<b>Uniprot ID</b>	<a href="#">P09104</a>
<b>NCBI</b>	<a href="#">NM_001975</a>
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...