


Product Datasheet

Villin (VIL1) antibody (orb1311338)

Description	Mouse monoclonal antibody to Villin (VIL1)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ChIP, IHC
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-188 of human VIL1 (Villin) (NP_009058) produced in E.coli.
Target	VIL1
Note	For research use only
Application notes	For reconstitution, we recommend adding 100uL distilled water to a final antibody concentration of about 1 mg/mL. To use this carrier-free antibody for conjugation experiment, we strongly recommend performing another round of desalting process.
Isotype	IgG1
Clonality	Monoclonal
MW	92.5 kDa
Uniprot ID	P09327
NCBI	NM_007127
Dilution Range	IHC 1:100~1000
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...