




Product Datasheet

CD44 antibody (orb1311576)

Description	Mouse monoclonal antibody to CD44
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Full length human recombinant protein of human CD44(NP_000601) produced in HEK293T cell.
Target	CD44
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P16070
NCBI	NM_000610
Dilution Range	IHC 1:500
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of paraffin...



HEK293T cells were transfected with the

...