


## Product Datasheet

### BubR1 (BUB1B) antibody (orb1311661)

<b>Description</b>	Mouse monoclonal antibody to BubR1 (BUB1B)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Monkey
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IF, IHC, WB
<b>Immunogen</b>	Protein expressed in 293T cell transfected with human BUB1B expression vector
<b>Target</b>	BUB1B
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	119.5 kDa
<b>Uniprot ID</b>	<a href="#">O60566</a>
<b>NCBI</b>	<a href="#">NM_001211</a>
<b>Dilution Range</b>	WB 1:2000, IHC 1:50, IF 1:100, Flow 1:100, IF 1:100
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of extracts (35ug)...