


Product Datasheet

SARS-CoV-2 N Protein antibody (orb1311881)

Description	Rabbit polyclonal antibody to SARS-CoV-2 N Protein
Species/Host	Rabbit
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Full length recombinant SARS-CoV-2 N protein produced in Expi293F.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Expiration Date	12 months from date of receipt.



Western blot analysis of full length rec...