


## Product Datasheet

### SARS-CoV-2 N Protein antibody (orb1311883)

<b>Description</b>	Rabbit polyclonal antibody to SARS-CoV-2 N Protein
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Virus
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, ICC, IHC, WB
<b>Immunogen</b>	Full length recombinant SARS-CoV-2 N protein produced in E.coli.
<b>Target</b>	N Protein
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	47kDa
<b>Dilution Range</b>	ELISA 1:5000-1:10000, WB 1:2000, IHC 1:500-1:1000, ICC 1:500-1:1000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of full length rec...