


Product Datasheet

B Raf (BRAF) antibody (orb1311934)

Description	Rabbit polyclonal antibody to B Raf (BRAF)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human BRAF
Target	BRAF
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	84.3 kDa
Uniprot ID	P15056
NCBI	NM_004333
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...