


Product Datasheet

MRP2 (ABCC2) antibody (orb1311941)

Description	Rabbit polyclonal antibody to MRP2 (ABCC2)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human ABCC2
Target	ABCC2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	174 kDa
Uniprot ID	Q92887
NCBI	NM_000392
Dilution Range	WB: 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...