


Product Datasheet

Cyclin A2 (CCNA2) antibody (orb1311953)

Description	Rabbit polyclonal antibody to Cyclin A2 (CCNA2)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human CCNA2
Target	CCNA2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	48.4 kDa
Uniprot ID	P20248
NCBI	NM_001237
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...