


Product Datasheet

HSPC142 (BABAM1) antibody (orb1311958)

Description	Rabbit polyclonal antibody to HSPC142 (BABAM1)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human BABAM1
Target	BABAM1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	36.4 kDa
Uniprot ID	Q9NWW8
NCBI	NM_001033549
Dilution Range	WB: 1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...