


Product Datasheet

ASK1 (MAP3K5) antibody (orb1311988)

Description	Rabbit polyclonal antibody to ASK1 (MAP3K5)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Recombinant protein of human MAP3K5
Target	MAP3K5
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	154.4 kDa
Uniprot ID	Q99683
NCBI	NM_005923
Dilution Range	WB: 1:500~1:1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...