


Product Datasheet

Fatty Acid Binding Protein 5 (FABP5) antibody (orb1311999)

Description	Rabbit polyclonal antibody to Fatty Acid Binding Protein 5 (FABP5)
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Synthetic peptide conjugated to KLH derived from C terminal of Human FABP5.
Target	FABP5
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	15 kDa
Uniprot ID	Q01469
NCBI	NM_001444
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...