

## Product Datasheet

### IL17 (IL17A) antibody (orb1312018)

|                            |   |
|----------------------------|---|
| <b>Description</b>         | Mouse monoclonal antibody to IL17 (IL17A)   |
| <b>Species/Host</b>        | Mouse   |
| <b>Reactivity</b>          | Human   |
| <b>Conjugation</b>         | Unconjugated  |
| <b>Tested Applications</b> | ELISA   |
| <b>Immunogen</b>           | Human recombinant protein fragment corresponding to amino acids 24-155 of human IL-17A (NP_002181) produced in E.coli.                  |
| <b>Target</b>              | IL17A   |
| <b>Storage</b>             | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| <b>Note</b>                | For research use only   |
| <b>Isotype</b>             | IgG1  |
| <b>Clonality</b>           | Monoclonal  |
| <b>MW</b>                  | 17.5 kDa  |
| <b>Uniprot ID</b>          | <a href="#">Q16552</a>  |
| <b>NCBI</b>                | <a href="#">NM_002190</a>   |
| <b>Dilution Range</b>      | ELISA 1:5000-10000  |
| <b>Expiration Date</b>     | 12 months from date of receipt.   |