

Product Datasheet

IL17 (IL17A) antibody (orb1312020)

Description	Mouse monoclonal antibody to IL17 (IL17A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Human recombinant protein fragment corresponding to amino acids 24-155 of human IL-17A (NP_002181) produced in E.coli.
Target	IL17A
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	17.5 kDa
Uniprot ID	Q16552
NCBI	NM_002190
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.