


Product Datasheet

SARS-CoV-2 RBD Neutralizing antibody (orb1312078)

Description	Mouse monoclonal antibody to SARS-CoV-2 RBD Neutralizing
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, NeA
Immunogen	Recombinant fusion protein of SARS-CoV-2 Spike RBD protein and virion surface domain of membrane protein (M-protein) produced in HEK293
Target	SARS-CoV-2 Spike Protein
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	34kDa
Dilution Range	ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.



ELISA using SARS-CoV-2 RBD Neutralizing ...