


Product Datasheet

SARS-CoV-2 S Protein antibody (orb1312120)

Description	Mouse monoclonal antibody to SARS-CoV-2 S Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant SARS-CoV-2 Spike-RBD protein produced in E.coli.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB 1:500-1000, ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant SAR...