


Product Datasheet

SARS-CoV-2 N Protein antibody (orb1312143)

Description	Mouse monoclonal antibody to SARS-CoV-2 N Protein
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant full length SARS-CoV-2 N protein produced in E.coli.
Target	N Protein
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Dilution Range	WB 1:500-1000, ELISA 1:5000-10000
Expiration Date	12 months from date of receipt.



Western blot analysis of recombinant full length SARS-CoV-2 N protein produced in E.coli.