


Product Datasheet

Peroxiredoxin 6 (PRDX6) antibody (orb1312486)

Description	Mouse monoclonal antibody to Peroxiredoxin 6 (PRDX6)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human Prdx6 (NP_004896) produced in E.coli.
Target	PRDX6
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	25 kDa
Uniprot ID	P30041
NCBI	NM_004905
Dilution Range	WB 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...