


## Product Datasheet

### Peroxiredoxin 6 (PRDX6) antibody (orb1312487)

<b>Description</b>	Mouse monoclonal antibody to Peroxiredoxin 6 (PRDX6)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Full length human recombinant protein of human Prdx6 (NP_004896) produced in E.coli.
<b>Target</b>	PRDX6
<b>Note</b>	For research use only
<b>Isotype</b>	IgG2b
<b>Clonality</b>	Monoclonal
<b>MW</b>	25 kDa
<b>Uniprot ID</b>	<a href="#">P30041</a>
<b>NCBI</b>	<a href="#">NM_004905</a>
<b>Dilution Range</b>	WB 1:500-1000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...