


Product Datasheet

HSP27 (HSPB1) antibody (orb1312519)

Description	Mouse monoclonal antibody to HSP27 (HSPB1)
Species/Host	Mouse
Reactivity	Human, Monkey
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human HSPB1 (NP_001531) produced in E.coli.
Target	HSPB1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2b
Clonality	Monoclonal
MW	22.8 kDa
Uniprot ID	P04792
NCBI	NM_001540
Dilution Range	WB 1:500-1000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...