


Product Datasheet

Nucleostemin (GNL3) antibody (orb1312561)

Description	Mouse monoclonal antibody to Nucleostemin (GNL3)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human GNL3 (NP_055181) produced in E.coli.
Target	GNL3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	61.8 kDa
Uniprot ID	Q9BVP2
NCBI	NM_014366
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...