

Product Datasheet

Calreticulin (CALR) antibody (orb1312605)

Description Mouse monoclonal antibody to Calreticulin (CALR)

Species/Host Mouse

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications WB

Immunogen Full length human recombinant protein of human CALR (NP_004334) produced in HEK293T cell.

Target CALR

Storage Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.

Note For research use only

Isotype IgG1

Clonality Monoclonal


MW 48.14 kDa

Uniprot ID [P27797](#)

NCBI [NM_004343](#)

Dilution Range WB 1:500

Expiration Date 12 months from date of receipt.



Western blot analysis of HEK293T cells w...