


Product Datasheet

PTER antibody (orb1312648)

Description	Mouse monoclonal antibody to PTER
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-242 of human PTER (NP_001001484) produced in E.coli.
Target	PTER
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38.8 kDa
Uniprot ID	Q96BW5
NCBI	NM_001001484
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...