


Product Datasheet

Diazepam Binding Inhibitor (DBI) antibody (orb1312692)

Description	Mouse monoclonal antibody to Diazepam Binding Inhibitor (DBI)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-104 of human DBI (NP_065438) produced in E.coli.
Target	DBI
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	11.6 kDa
Uniprot ID	P07108
NCBI	NM_020548
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...