


Product Datasheet

Stefin B (CSTB) antibody (orb1312710)

Description	Mouse monoclonal antibody to Stefin B (CSTB)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 1-98 of human CSTB (NP_000091) produced in E.coli.
Target	CSTB
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	11 kDa
Uniprot ID	P04080
NCBI	NM_000100
Dilution Range	WB 1:500, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...