


Product Datasheet

PCDH15 antibody (orb1312779)

Description	Mouse monoclonal antibody to PCDH15
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 160-400 of human PCDH15 (NP_001136235) produced in E.coli.
Target	PCDH15
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	217.06 kDa
Uniprot ID	Q96QU1
NCBI	NM_001142763
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...