


Product Datasheet

SSX4 (SSX4B) antibody (orb1312892)

Description	Mouse monoclonal antibody to SSX4 (SSX4B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SSX4B (NP_001030004) produced in E.coli.
Target	SSX4B
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O60224
NCBI	NM_001034832
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...