


## Product Datasheet

### MRP2 (ABCC2) antibody (orb1312946)

<b>Description</b>	Mouse monoclonal antibody to MRP2 (ABCC2)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 1256-1545 of human ABCC2 (NP_000383) produced in E.coli.
<b>Target</b>	ABCC2
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">Q92887</a>
<b>NCBI</b>	<a href="#">NM_000392</a>
<b>Dilution Range</b>	WB 1:500, IHC 1:5000
<b>Expiration Date</b>	12 months from date of receipt.



Western blot analysis of HEK293T cells w...