


Product Datasheet

Prothrombin (F2) antibody (orb1312952)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to Prothrombin (F2) |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 44-622 of human F2 (NP_000497) produced in HEK293T cell. |
| Target | F2 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG2a |
| Clonality | Monoclonal |
| MW | 70.01 kDa |
| Uniprot ID | P00734 |
| NCBI | NM_000506 |
| Dilution Range | WB 1:500 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...