


Product Datasheet

gamma Catenin (JUP) antibody (orb1312987)

Description	Mouse monoclonal antibody to gamma Catenin (JUP)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 500-745 of human JUP (NP_002221) produced in E.coli.
Target	JUP
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	81.6 kDa
Uniprot ID	P14923
NCBI	NM_002230
Dilution Range	WB 1:500, IHC 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...