


Product Datasheet

Genethonin 1 (STBD1) antibody (orb1313044)

Description	Mouse monoclonal antibody to Genethonin 1 (STBD1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 24-358 of human STBD1 (NP_003934) produced in E.coli.
Target	STBD1
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	38.8 kDa
Uniprot ID	O95210
NCBI	NM_003943
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...