

## Product Datasheet

### Genethonin 1 (STBD1) antibody (orb1313045)

<b>Description</b>	Mouse monoclonal antibody to Genethonin 1 (STBD1)
<b>Species/Host</b>	Mouse
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	WB
<b>Immunogen</b>	Human recombinant protein fragment corresponding to amino acids 24-358 of human STBD1 (NP_003934) produced in E.coli.
<b>Target</b>	STBD1
<b>Note</b>	For research use only
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>MW</b>	38.8 kDa
<b>Uniprot ID</b>	<a href="#">O95210</a>
<b>NCBI</b>	<a href="#">NM_003943</a>
<b>Dilution Range</b>	WB 1:500
<b>Expiration Date</b>	12 months from date of receipt.

