


Product Datasheet

Angiotensin Converting Enzyme 1 (ACE) antibody (orb1313086)

Description	Mouse monoclonal antibody to Angiotensin Converting Enzyme 1 (ACE)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 642-1230 of human ACE (NP_000780) produced in E.coli.
Target	ACE
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	P12821
NCBI	NM_000789
Dilution Range	WB 1:500, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...