


Product Datasheet

Arp3 (ACTR3) antibody (orb1313095)

Description	Mouse monoclonal antibody to Arp3 (ACTR3)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 2-418 of human ACTR3 (NP_005712) produced in E.coli.
Target	ACTR3
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47.2 kDa
Uniprot ID	P61158
NCBI	NM_005721
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of extracts (35ug)...