


Product Datasheet

BRUNOL6 (CELF6) antibody (orb1313208)

Description	Mouse monoclonal antibody to BRUNOL6 (CELF6)
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CELF6 (NP_443072) produced in E.coli.
Target	CELF6
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	50.3 kDa
Uniprot ID	Q96J87
NCBI	NM_052840
Dilution Range	WB 1:250~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...