


Product Datasheet

UBE1C (UBA3) antibody (orb1313235)

| | |
|----------------------------|---|
| Description | Mouse monoclonal antibody to UBE1C (UBA3) |
| Species/Host | Mouse |
| Reactivity | Canine, Human, Monkey, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 1-463 of human UBA3 (NP_003959) produced in E.coli. |
| Target | UBA3 |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 51.7 kDa |
| Uniprot ID | Q8TBC4 |
| NCBI | NM_003968 |
| Dilution Range | WB 1:500 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...