


Product Datasheet

Ovary specific acidic protein (MGARP) antibody (orb1313286)

Description	Mouse monoclonal antibody to Ovary specific acidic protein (MGARP)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human MGARP (NP_116012) produced in E.coli.
Target	MGARP
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q8TDB4
NCBI	NM_032623
Dilution Range	WB 1:2000, IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...