

Product Datasheet

SPANXB1 antibody (orb1313530)

Description	Mouse monoclonal antibody to SPANXB1
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human SPANXB1 (NP_115850) produced in E.coli.
Target	SPANXB1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	Q9NS25
NCBI	NM_032461
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.

