


Product Datasheet

p107 (RBL1) antibody (orb1313561)

| | |
|----------------------------|--|
| Description | Mouse monoclonal antibody to p107 (RBL1) |
| Species/Host | Mouse |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Human recombinant protein fragment corresponding to amino acids 739-1014 of human RBL1 (NP_899662) produced in E.coli. |
| Target | RBL1 |
| Note | For research use only |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| MW | 114.5 kDa |
| Uniprot ID | P28749 |
| NCBI | NM_183404 |
| Dilution Range | WB 1:500 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of HEK293T cells w...