


Product Datasheet

KRT84 antibody (orb1313659)

Description	Mouse monoclonal antibody to KRT84
Species/Host	Mouse
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human KRT84 (NP_149034) produced in E.coli.
Target	KRT84
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	64.7 kDa
Uniprot ID	Q9NSB2
NCBI	NM_033045
Dilution Range	IHC 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...