


Product Datasheet

TCP1 delta (CCT4) antibody (orb1313875)

Description	Mouse monoclonal antibody to TCP1 delta (CCT4)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human CCT4 (NP_006421) produced in E.coli.
Target	CCT4
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	57.7 kDa
Uniprot ID	P50991
NCBI	NM_006430
Dilution Range	WB 1:500~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...