


Product Datasheet

C3orf43 (SMCO1) antibody (orb1315461)

Description	Mouse monoclonal antibody to C3orf43 (SMCO1)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human C3orf43 (NP_001071125) produced in E.coli.
Target	SMCO1
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	24.4 kDa
Uniprot ID	Q147U7
NCBI	NM_001077657
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...