


Product Datasheet

TRBP (TARBP2) antibody (orb1315586)

Description	Mouse monoclonal antibody to TRBP (TARBP2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Full length human recombinant protein of human TARBP2 (NP_004169) produced in E.coli.
Target	TARBP2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	36.7 kDa
Uniprot ID	Q15633
NCBI	NM_004178
Dilution Range	WB 1:2000, IHC 1:150
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...