


Product Datasheet

Troponin I fast skeletal muscle (TNNI2) antibody (orb1315794)

Description	Mouse monoclonal antibody to Troponin I fast skeletal muscle (TNNI2)
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Full length human recombinant protein of human TNNI2 (NP_003273) produced in E.coli.
Target	TNNI2
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	21.2 kDa
Uniprot ID	P48788
NCBI	NM_003282
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...