


Product Datasheet

BAFF (TNFSF13B) antibody (orb1315807)

Description	Mouse monoclonal antibody to BAFF (TNFSF13B)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 134-285 of human TNFSF13B(NP_006564) produced in E.coli.
Target	TNFSF13B
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
MW	31 kDa
Uniprot ID	Q9Y275
NCBI	NM_006573
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...