


Product Datasheet

VPS13B antibody (orb1315808)

Description	Mouse monoclonal antibody to VPS13B
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 91-346 of human VPS13B(NP_858047) produced in E.coli.
Target	VPS13B
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
MW	47 kDa
Uniprot ID	Q7Z7G8
NCBI	NM_181661
Dilution Range	WB 1:2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...