


Product Datasheet

DR4 (TNFRSF10A) antibody (orb1316002)

Description	Mouse monoclonal antibody to DR4 (TNFRSF10A)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 263-468 of human TNFRSF10A(NP_003835) produced in E.coli.
Target	TNFRSF10A
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Uniprot ID	O00220
NCBI	NM_003844
Dilution Range	WB 1:200~2000
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...