


Product Datasheet

IL23 Receptor (IL23R) antibody (orb1316052)

Description	Mouse monoclonal antibody to IL23 Receptor (IL23R)
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	Human recombinant protein fragment corresponding to amino acids 57-355 of human IL23R(NP_653302) produced in E.coli.
Target	IL23R
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q5VWK5
NCBI	NM_144701
Dilution Range	WB 1:500
Expiration Date	12 months from date of receipt.



Western blot analysis of HEK293T cells w...